




# FURNACE ZE



**LABORATORY FURNACE**

*- Laboratory Sector -*



## LABORATORY FURNACE

The ZE electric muffle furnace, more compact than the Z series, is used in laboratories, specifically research, educational, industrial and analytical laboratories.

The ZE oven is heated on four sides by a muffle complete with heating elements to guarantee maximum temperature uniformity. It is also finely painted and has various safety systems: microswitch and magnetothermal to protect the operator and the electrical system itself.

The furnace has a digital thermoregulator, a static unit and, on request, a temperature programmer, daily or weekly timer, from which it is possible to set automatic on/off.





## DATA SHEET

### TECHNICAL DATA

TYPE	SINGLE-PHASE VOLTAGE	FREQUENCY	WATTAGE	MAX TEMPERATURE
ZE	v.110 – v.230	Hz. 50-60	2.3 kw	1100 C°

### OVEN CHAMBER DIMENSIONS

TYPE	DEPTH	WIDTH	HEIGHT
ZE	250 mm	145 mm	100 mm

### DIMENSIONS

TYPE	DEPTH	WIDTH	HEIGHT	TYPE
ZE	580 mm	400 mm	540 mm	40 kg